

ESfEQA REPORT INTERPRETATION

EXPLANATION MANUAL FOR SURVEY REPORTS
AND STATISTICAL EVALUATIONS



ESfEQA PROGRAMS



BIOCHEMISTRY

Blood Gas and Electrolytes, Cardiac Marker, Cerebrospinal Fluid, Clinical Chemistry, Coagulation, Drugs of Abuse, Ethanol, Fecal Occult Blood, Glycated Hemoglobin, Qualitative Urine Analysis, Therapeutic Drugs, Urine Chemistry, Urine Sediments

IMMUNOLOGY

hCG, Hormones, Procalcitonin, Specific Proteins, Tumor Marker

MICROBIOLOGY

Adenovirus, Aspergillus Fumigatus, Bacteriology (Pathogen Identification and Antimicrobial Susceptibility Testing according to EUCAST and CLSI), Borrelia, Brucella, Chagas, Chlamydomphila Pneumoniae, Chlamydomphila Trachomatis, Chikungunya Virus, Coxsackie-virus, Dengue Virus, ECHO Virus, Enterovirus, Epstein-Barr Virus, Hepatitis A Virus, Hepatitis B Virus, Hepatitis E Virus, HTLV I/II, Infectious Disease (HIV, HCV, HBV), Influenza A Virus, Influenza B Virus, Leptospira, Malaria Microscopy (Parasite/Stage Identification), Measles, Mosquito Transmitted Diseases (Chikungunya, West-Nile, Zika, Dengue), Parainfluenza Virus, Parvovirus B 19, Respiratory Syncytial Virus, SARS-CoV-2, Syphilis, ToRCH, Varicella Zoster Virus, West-Nile-Fever Virus, Zika Virus

HEMATOLOGY

Blood Grouping, Immunhematology, Erythrocyte Sedimentation Rate, Erythrocyte Sedimentation Rate for Alifax and Alcor devices, Hemogram, Hemogram including 3-part Diff.

EXTERNAL QUALITY ASSESSMENT BY ESfEQA

Quality Standards in Medical Laboratories

A high quality standard is essential for every medical laboratory since test results are the basis for medical decisions and have an important impact on the well-being and treatment of patients.

There are different approaches to maintain and improve quality in medical laboratories. ESfEQA -European Society for External Quality Assessment - supports laboratories to assess the quality of their analytical results and ultimately improve their performance by providing well-designed external quality assessment programs.

ESfEQA offers a wide range of External Quality Assessment Schemes

ESfEQA was founded in Heidelberg/Germany in 2013 and is accredited according to the international standard ISO 17043:2010 by the German national accreditation body DAkkS. Since its foundation, ESfEQA has expanded its program portfolio and at the same time the number of laboratories participating in ESfEQA's external quality assessment schemes. Currently, ESfEQA offers more than 70 quantitative and qualitative EQA programs worldwide in the areas of biochemistry, immunology, microbiology and hematology.


Participants submit their results to the ESfEQA database via the secured TEQA web application. ESfEQA evaluates the participants' results according to ISO 13528 and provides the reports to the participants via TEQA.

External Quality Assessment Reports

This manual is intended to explain the survey reports and statistical evaluations to the participants. ESfEQA reports are designed to provide participants with detailed information on their laboratory performance in order to assess their level of quality and to show how the the quality of the laboratory performance can be improved if quality standards are not fully met. ESfEQA reports provide a comparison of laboratory results with other laboratories that use the same analytical system, in particular the same analytical device or reagent, but also provide detailed information for comparison to other analytical devices. We consider the comparison to other analytical devices, methods etc. to be very important since the medical interpretation of laboratory results is usually not associated to any particular instrument, method or reagent brand.

In addition to detailed information, ESfEQA reports contain a summary that allows participants to capture the laboratory performance at a glance.

We are looking forward to support you with our enthusiasm for quality in medical laboratories in maintaining and improving the quality standard in your laboratory.



Summary Page

The summary page lists all analytes that were determined by the participant and that were statistically evaluated. For each analyte, the participant result, z-score and bias are displayed.

Clinical Chemistry 4										
1 4 (28/10/2019) DISPATCH					2 11/11/2019 09:51:22 / 1 PROCESSING DATE / VERSION					
Report Summary										
3 sample: C06019_04 8 10 11 12 13										
Parameter	5 Evaluation	Target value	SD	U	Number of Labs results		Lab result	Z-Score	BIAS [%]	
Cholesterol mg/dL	GENERAL	195 (A)	8.45 (S)	0.51	430	430		-1.3	-5.64	
	INSTRUMENT GROUP	196 (A)	8.1 (S)	1.66	41	41	184	-1.5	-6.36	
	REAGENT	202 (A)	8.1 (S)	2.93	14	14		-2.1	-9.18	
Creatinine mg/dL	GENERAL	3.802 (A)	0.24 (S)	0.01	450	450		-1.4	-9	
	INSTRUMENT GROUP - METHOD	3.998 (A)	0.25 (S)	0.05	35	35	346	-2.1	-13.46	
	METHOD	3.787 (A)	0.24 (S)	0.02	304	304		-1.4	-8.63	
Glucose mg/dL	GENERAL	258.6 (A)	12.93 (S)	0.78	425	425		1.7	8.7	
	INSTRUMENT GROUP	258.9 (A)	12.94 (S)	2.59	39	39	281	1.7	8.57	
	REAGENT	268.5 (A)	13.43 (S)	4.33	15	15		0.9	4.69	
ALP Alkaline Phosphatase U/L	GENERAL	224.3 (A)	16.45 (S)	3.3	394	394		5.7	41.77	
	INSTRUMENT GROUP	269.9 (A)	19.79 (S)	9.53	38	38	318	2.4	17.82	
	METHOD	280.3 (A)	20.56 (S)	13.76	31	31		1.8	13.45	
Total Protein g/dL	GENERAL	7.385 (A)	0.25 (S)	0.02	371	371		4.1	13.74	
	INSTRUMENT GROUP	7.504 (A)	0.25 (S)	0.05	36	36	8.4	3.6	11.94	
	REAGENT	8.007 (A)	0.27 (S)	0.1	11	11		1.5	4.91	
Triglycerides mg/dL	GENERAL	212.5 (A)	11.33 (S)	0.69	426	426		-0.2	-1.18	
	INSTRUMENT GROUP	205.7 (A)	10.97 (S)	2.14	41	41	210	0.4	2.09	
	REAGENT	214.6 (A)	11.44 (S)	3.69	15	15		-0.4	-2.14	
Urea mg/dL urea N	GENERAL	67.05 (A)	4.47 (S)	0.27	432	432		2	13.04	
	INSTRUMENT GROUP	68.59 (A)	4.57 (S)	0.88	42	42	75.8	1.6	10.5	
	REAGENT	64.86 (A)	4.32 (S)	1.27	18	18		2.5	16.85	
Uric Acid mg/dL	GENERAL	9.445 (A)	0.41 (S)	0.03	408	408		-0.6	-2.59	
	INSTRUMENT GROUP	9.832 (A)	0.43 (S)	0.08	40	40	9.2	-1.5	-6.43	
	REAGENT	10.24 (A)	0.44 (S)	0.14	16	16		-2.3	-10.16	
Bilirubin Direct mg/dL	GENERAL	2.236 (A)	0.16 (S)	0.01	376	376		-1.9	-13.69	
	INSTRUMENT GROUP	2.343 (A)	0.17 (S)	0.04	36	36	1.93	-2.4	-17.63	
	REAGENT	2.216 (A)	0.16 (S)	0.05	14	14		-1.8	-12.91	
Bilirubin Total mg/dL	GENERAL	5.018 (A)	0.37 (S)	0.02	417	417		2	14.59	
	INSTRUMENT GROUP	5.077 (A)	0.37 (S)	0.08	38	38	5.75	1.8	13.26	
	REAGENT	5.24 (A)	0.38 (S)	0.13	13	13		1.3	9.73	

* (A) AlgoA, (G) Expert Group, (R) Reference value, (C) Certified Reference, (S) Set by coordinator, (E) Experiment, (I) Interquartile, (N) Non robust sd

IMPORTANT: The z-score is the most important indicator for the performance of the laboratory.

- A z-score of 0 represents a participant result that is identical to the target value.
- A z-score of +3 indicates a participant results that is at the upper limit of the permissible range
- A z-score of -3 indicates a participant results that is at the lower limit of the permissible range
- Certificates of successful participation are issued if the z-score of at least one of the calculated groups of an analyte is within the range of ± 3 for all samples of the particular survey.

Header

- 1 Dispatch number: the dispatches are counted continuously. The deadline for result submission is displayed in brackets.
- 2 The processing data and version indicate the date when the statistical evaluation for this survey was performed and which version of the statistical evaluation is provided. Usually, there is only a single version (version number "1") of the statistical evaluation.

Report Summary


- 3 Sample identifier
- 4 Units: units indicated in the report reflect the units that were used by the participant for result submission.
- 5 Evaluated groups: in this example for Cholesterol, the evaluation of the general group, the instrument group, and the reagent group is provided. The general group includes results from all participants that submitted results for this analyte. The instrument (reagent) group summarizes results from all laboratories that have used the same instrument (reagent) as the participant of this report. ESfEQA predefines appropriate evaluation groups for each analyte.
- 6 Target values of the evaluated groups. Target values can be group specific, e.g. the target value for the general group might be different from the instrument group and the target value for different instrument groups might be instrument specific.
- 7 Statistical basis for the determination of target values:
 - (A) = calculated as robust mean (Algorithm A) according to ISO 13528
 - (R) = reference value determined by an expert laboratory
 - (C) = certified reference value of certified material
 - (D) = calculated as robust mean of a group of expert laboratories
- 8 SD = Standard Deviation (Definition of SD see below)
- 9 Statistical basis how the SD was determined:
 - (S) = set by coordinator:
 - SD = (permissible deviation [%] / 100) / 3 · X if the relative deviation applies and
 - SD = (permissible deviation / 3) if the absolute deviation applies.
 - (E) = determined by experiment

(I) = determined as interquartile range
 (N) = determined as non-robust SD
 In most cases of ESfEQA surveys, the permissible deviation for each analyte is derived from its medical relevance and the reference interval for the analyte. For most quantitative analytes ESfEQA, in collaboration with its scientific advisors, has defined a relative permissible deviation (in %) and below a defined concentration threshold an absolute permissible deviation (expressed in the unit of the analyte, e.g. U/L, mmol/l).

- 10 U = Uncertainty of the assigned value, calculated as:
 $U = 1,25 \cdot s^* / \sqrt{p}$ (s* as robust standard deviation calculated by the Algorithm A and p as number of results)
 An uncertainty above $0,3 \cdot SD$ is not negligible and might result in an inappropriate calculation of the z-score.
- 11 Number of laboratories that submitted results for this analyte and total number of results submitted by these laboratories
- 12 Z-score, calculated as:
 $ZS = (x-X)/SD$ (x as participant value, X as target value)
 For results that were included in several groups (e.g. I, M group and general group) there are several target values X and thus several z-scores.
- 13 BIAS [%] = Deviation of the participants result from the target value in [%]
- 14 this symbol indicates a z-score above 2,0 or below -2,0 (in this case -2,1 in the instrument-method group) and shall be considered as a warning signal. It is displayed in the line of an analyte if any of the z-cores of the evaluated groups is outside of this range.
- 15 this symbol indicates a z-score above 3,0 or below -3,0 (in this case 5,7 in the general group) and shall be considered as an action signal. It alerts the participant about instrument-, method and reagent-specific target values. It is displayed in the line for an analyte if any of the z-cores of the evaluated groups is outside of this range.

Analyte Specific Information (quantitative surveys)

In this section of the report, a Frequency Histogram, Statistical Indices, Youden Diagram and Statistical Comparison Data are provided.



Header

The HEADER contains the laboratory code and laboratory name, the name of the analyte, the default unit, the dispatch of the EQA survey, the laboratory configuration, the date of the statistical evaluation and the version of data processing.

- Unit used by the participant for result submission.
- Lab configuration: Analyte-specific configuration of the participating laboratory. The evaluated participant groups that are displayed in the charts for this analyte are those defined by the participant's configuration. In this case, these groups are:
 - General (including the results of all instrument and reagent groups)
 - Instrument group: Beckmann Coulter AU
Several comparable instruments can be combined into one group. In this case, the 'Beckmann Coulter AU instrument' was evaluated in the broad instrument group 'Beckmann Coulter AU'.
 - Reagent: Beckmann Coulter System-Reagent

Frequency Histograms

The FREQUENCY HISTOGRAMS indicate the result distributions of the various groups (e.g. general and instrument/method) that were evaluated.

- Groups in which the result of the participant was evaluated, in this case:
 - General
 - Instrument group: Beckmann Coulter AU
 - Reagent: Beckmann Coulter System-Reagent
- Number of results submitted for the selected analyte (in this example, 'Albumin') in the reagent group "Beckmann Coulter System-Reagent" that are between 2,509 g/dL and 2,609 g/dL.
- Number of results submitted for the selected analyte (in this example, 'Albumin') in the instrument group "Beckmann Coulter AU" that are between 2,609 g/dL and 2,709 g/dL.
- Number of results submitted for the selected analyte (in this example, 'Albumin') in the general group that are between 2,809 g/dL and 2,909 g/dL.
- The check mark indicates the location of the value submitted by the participant.

Statistical Indices

The STATISTICAL INDICES provide information about the target value and how it was calculated (e.g. as participant consensus calculated by the Algorithm A), the SD (Standard deviation), the uncertainty of the assigned value, the number of labs and results that were considered for the statistical evaluation. The last three columns display the lab result, the z-score (ZS) and the BIAS (deviation of the lab result to the target value in percent).

- Groups in which the result of the participant was evaluated, in this case:
 - General
 - Instrument group: Beckmann Coulter AU
 - Reagent: Beckmann Coulter System-Reagent
- CV = Coefficient of Variation, calculated as:
 $CV = SD/X \cdot 100$

Youden Diagram

(only displayed for EQAs with two samples per survey)

- Blue dot: Display of the participant's z-score for the selected analyte (in this example, 'Albumin') determined in sample CC4_2019_04_a compared to the z-score for the selected analyte determined in sample CC4_2019_04_b.
- Yellow dots represent the z-scores of all other participants (of the evaluated group) for the selected analyte (in this example, 'Albumin') determined in sample CC4_2019_04_a compared to the z-score for the selected analyte determined in sample CC4_2019_04_b.


Analyte Specific Information (quantitative surveys)

Statistical
Comparison

Statistical Comparison							
Evaluation	Target value	SD		Uncertainty *	CV	Number of	
	AlgoA	set by coordinator				Labs	results
1 general	2.859	0.2	-	0.01	7	365	365
INSTRUMENT GROUP							
2 Roche Cobas	2.894	0.2	-	0.03	6.91	95	95
3 Beckmann Coulter AU	2.697	0.2	-	0.03	7.42	61	61
4 Mindray BS Series	2.939	0.2	-	0.04	6.81	31	31
5 ORTHO Vitros	9.886	0.66	-	0.24	6.67	12	12
6 BIOSYSTEMS	2.822	0.2	-	0.08	7.09	9	9
7 ABBOTT Architect	2.823	0.2	-	0.05	7.08	25	25
8 HUMAN HumaLyzer Photometer	2.883	0.2	-	0.06	6.94	19	19
9 HUMAN HumaStar Analyzers	2.795	0.2	-	0.05	7.16	24	24
10 ELITECHGroup Selectra Pro	3.122	0.21	-	0.12	6.67	5	5
11 ERBA Mannheim	2.993	0.2	-	0.09	6.68	8	8
12 Roche Cobas c111	3.042	0.2	-	0.09	6.67	8	8
REAGENT							
13 ROCHE System-Reagent	2.906	0.2	-	0.03	6.88	96	96
14 BECKMAN COULTER System Reagent	2.692	0.2	-	0.03	7.43	52	52
15 HUMAN Reagent	2.794	0.2	-	0.04	7.16	47	47
16 ABBOTT System-Reagent	2.826	0.2	-	0.05	7.08	21	21
17 ORTHO System-Reagent/Slide	9.886	0.66	-	0.24	6.67	12	12
18 SPINREACT Clinical Chemistry	3.119	0.21	-	0.08	6.67	10	10
19 MINDRAY System-Reagent	2.912	0.2	-	0.08	6.87	9	9
20 ELITECH Group Reagents	3.104	0.21	-	0.12	6.67	5	5
21 DIASYS System Reagent	2.69	0.2	-	0.11	7.43	5	5
22 BIOSYSTEMS System Reagent	2.836	0.2	-	0.11	7.05	5	5
23 ERBA MANNHEIM System-Reagent	3.106	0.21	-	0.12	6.67	5	5
24 SIEMENS System-Reagent	2.75	0.2	-	0.11	7.27	5	5

* Uncertainty of the assigned value

Shewhart
Chart



Statistical Comparison

The table STATISTICAL COMPARISON indicates the target value, the SD, the uncertainty (of the calculated target value), the CV and the number of participating laboratories and results for each evaluated group.

Shewhart Chart

For each laboratory, the SHEWHART CHART displays the z-scores of up to 12 previous surveys for the selected analyte of the EQA scheme. The sample number, the deadline of result entry, the sample identifier, the target value, the z-scores and the number of labs/results are displayed in the table. Numbers are provided for all groups (e.g. general, instrument/method) that were evaluated.

- 12 Blue dots: Each blue dot in the SHEWHART Chart indicates the z-score of the participant for the selected analyte (in this example, 'Albumin') and the specific comparison group (in this example, the general group) that was determined in the corresponding survey sample:
 - Dispatch 9
 - Sample CC4_2018_01_a
 - Target value of the general group: 3,462 g/dL
 - Z-score: -1,1
 - Number of labs: 102/number of results: 102

- 13 Blue triangles: Each blue triangle indicates the z-score of the participant for the selected analyte (in this example, 'Amylase') obtained in the second sample of a survey (only displayed for EQAs with two samples per survey):
 - Dispatch 9
 - Sample CC4_2018_01_b
 - Target value of the general group: 5,247 g/dL
 - Z-score: -0,4
 - Number of labs: 151/number of results: 151

Analyte Specific Information (qualitative surveys)

A table of results and a bar diagram illustrate the results for each analyte.




IMPORTANT: A certificate of successful participation for each individual analyte is issued if the lab results for all samples of this survey are identical to the expected result (true result)

Analyte Specific Information (qualitative surveys)

The TABLE indicates the results of all participants for each evaluated group (e.g. general and instrument/reagent)

- 1 True result = Expected result defined by ESFEQA, usually derived from the results of predistribution testings of the sample
- 2 Lab result = participants result
- 3 The tick indicates the participant's result
- 4 Total number of results submitted as "positive"
- 5 Total number of results submitted as "borderline"
- 6 Total number of results submitted as "negative"

Cycle Report



A table of results and a bar diagram illustrate the results for each analyte.

Header

The HEADER indicates the laboratory code and laboratory name, the analyte, the lab configuration, the selected period of time for evaluation and dates of the first and the last dispatch that were selected for the cycle report (range).

- 1 Evaluation period of the cycle report that has been selected by the participant in TEQA.
- 2 The range of dispatches indicates the date of the first and the last dispatch that included in the cycle evaluation.

Table

The TABLE in the cycle report displays the z-scores of all dispatches for an analyte in the evaluated groups over the selected time period in a color-coded manner.

- 3 Z-score for Fibrinogen, sample COA12_2019_08_a, deadline 02.09.2019, for the reagent group. The green color indicates that the absolute values of all z-scores in this evaluated group (reagent group) are below 2 ($|Zs| < 2$).
- 4 The line with the instrument group is not highlighted in any color-code since the Z-score for Fibrinogen is not consistently within the ± 2 range.
- 5 Red color indicates absolute z-scores above 3 ($|Zs| \geq 3$).
- 6 Yellow color indicates absolute z-scores between 2 and 3 ($2 \geq |Zs| < 3$).
- 7 Twelve dispatches were evaluated in this cycle report. This number is also used for the graph 'Laboratory performance comparison' below.
- 8 This number indicates the total number of results in the evaluated group.

Laboratory performance comparison

The LABORATORY PERFORMANCE COMPARISON charts reflect the performance of the laboratory compared to all other participants for each evaluated group over the dispatches of the selected time period. The z-score scale contains absolute values. Values exceeding ± 4 have been discarded.

- 9 First dispatch evaluated in this cycle report, sample COA12_2019_01, deadline 11.02.2019 (data indicated in the table above).

- 10 Lower border: average of the best 10% z-scores
- 11 Higher border: average of the worst 10% z-scores
- 12 Ascending arrow: indicates the absolute value of the z-score of the participant in the instrument group for the 5th evaluated dispatch (deadline 03.06.19, sample name COA12_2019_05_a), ascending arrow indicating a positive z-score.
- 13 Descending arrow: indicates the absolute value of the z-score of the participant in the instrument group for the 10th evaluated dispatch (deadline 28.10.19, sample name COA12_2019_10_a), descending arrow indicating a negative z-score.

ESFEQA GmbH

European Society for External
Quality Assessment
Siemensstr. 38
69123 Heidelberg
GERMANY
Phone +49(0)6221-894669-70
Fax +49(0)6221-894669-90
info@esfeqa.eu
www.esfeqa.eu

Authorized distributor

