

GLIDERITE STYLETS



DESIGNED TO WORK WITH GLIDESCOPE VIDEO LARYNGOSCOPES

The angle of GlideRite® Stylets complements the unique angle of GlideScope instruments to help facilitate quick placement of endotracheal tubes and help reduce patient trauma.

- Enhances ETT manoeuvrability for fast intubation
- Easy to high-level disinfect or sterilise in an autoclave
- Constructed of durable, reusable stainless steel
- Stylet and ETT insertion can be viewed in real time, helping ensure accurate tube placement

RIGID STYLET

- Designed for use in ET tubes 6.0 mm and larger
- Available in packs of 10
- Stylet rod length is 266 mm

DLT STYLET

- Designed for use in double lumen endotracheal tubes 6.0 mm and larger
- Available in packs of 10
- Stylet rod length is 394 mm

glidescope.com



verathon.com



Verathon Medical (Canada) ULC
2227 Douglas Road
Burnaby, BC V5C 5A9
Canada
Tel: 604 439 3009
Fax: 604 439 3039



Verathon Medical (Europe) B.V.
Willem Fenengastraat 13
1096 BL Amsterdam
The Netherlands
Tel: +31 (0) 20 210 30 91
Fax: +31 (0) 20 210 30 92



Verathon Inc.
20001 North Creek Parkway
Bothell, WA 98011
USA
Tel: 800 331 2313 (USA & Canada only)
Tel: 425 867 1348 Fax: 425 883 2896