



# ARIES<sup>®</sup>

ARIES<sup>®</sup> Systems | Welcome to the New Way to Work

## ARIES<sup>®</sup> Extraction Kit



## Welcome to the New Way to Work

---

**The ARIES<sup>®</sup> Extraction Kit enables laboratorians to take advantage of automation provided by ARIES<sup>®</sup> Systems. Designed for nucleic acid extraction and concentration, the ARIES<sup>®</sup> Extraction Kit offers:**

- Universal extraction of both DNA and RNA with a single cassette\*
- Internal sample processing controls to verify sample lysis, nucleic acid extraction, and proper system and reagent performance
- Walkaway processing, once the specimen is added to the sample chamber

\* Only validated with the sample types listed in Specifications.



**ARIES<sup>®</sup> Systems are crafted to increase laboratory efficiency, ensure result accuracy, and fit seamlessly into today's lean laboratory.**

## Specifications

Validated Specimen Types	Universal transport media (UTM), whole blood, stool, cerebrospinal fluid (CSF), bronchoalveolar lavage (BAL), plasma, and urine specimens. <i>Preprocessing steps may be required for some sample types.</i>
Input Sample Volume	200–400 µl
Required Instrumentation	ARIES® Systems

## Ordering Information

Product Name	Part Number
ARIES® Extraction Kit Includes: (24) Extraction Cassettes ARIES® Extraction Protocol File Package Insert Empty ARIES® PCR Tubes	50-10028
ARIES® Extraction Cassettes Reorder kit: 24 Cassettes	50-10026
ARIES® Two Module System Includes: Instrument System Operation Manual Two Magazines Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand	ARIES-M12V1-IVD
ARIES® M1 System Includes: Instrument System Operation Manual One Magazine Quick Guide Two Sample Prep Trays Power Cord Handheld Barcode Scanner and Stand	ARIES-M6V1-IVD

**Luminex**  
complexity simplified.

[support@luminexcorp.com](mailto:support@luminexcorp.com) or [orders@luminexcorp.com](mailto:orders@luminexcorp.com)

For In Vitro Diagnostic Use. For use with ARIES® System. The ARIES Extraction Kit has not been validated for use with any specific analytical test method. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at [support@luminexcorp.com](mailto:support@luminexcorp.com) to obtain the appropriate product information for your country of residence. ARIES Systems are class 1(I) laser products.

© 2016-2017 Luminex Corporation. All rights reserved. ARIES and Luminex are trademarks of Luminex Corporation, registered in the U.S. and other countries.

[www.luminexcorp.com](http://www.luminexcorp.com)

HEADQUARTERS

UNITED STATES

+1.512.219.8020

[info@luminexcorp.com](mailto:info@luminexcorp.com)

EUROPE

+31.73.800.1900

[europe@luminexcorp.com](mailto:europe@luminexcorp.com)

CANADA

+1.416.593.4323

[info@luminexcorp.com](mailto:info@luminexcorp.com)

CHINA

+86.21.802.31150

[info@luminexcorp.com](mailto:info@luminexcorp.com)

JAPAN

+81.3.5545.7440

[infojp@luminexcorp.com](mailto:infojp@luminexcorp.com)

SS19304.0317