

Luminex[®]

ARIES[®]

ARIES[®] Systems | Welcome to the New Way to Work

MultiCode[®] Ready Mix




ARIES[®] MultiCode[®] Ready Mixes provide the key PCR ingredients to run user-defined protocols on any ARIES[®] System. These mixes come in individual sealed tubes that users can simply add target and control primers to, and then attach directly to ARIES[®] Extraction Cassettes.

Features include:

- Lyophilized mixes that eliminate the need to prepare and thaw bulk PCR master mixes
- Ready to use reagents that help to improve efficiency and reduce costs
- Streamlined workflow that reduces preparation time



Workflow



Ordering Information*

Product Name	Part Number
ARIES® MultiCode® Ready Mix (Contains reverse transcriptase)	3696 (24 tubes)
ARIES® MultiCode® DNA Ready Mix	3697 (24 tubes)
ARIES® Two Module System (RUO)	ARIES-M12V1-RUO
ARIES® M1 System (RUO)	ARIES-M6V1-RUO
SYNCT™ Software	CN-SW34

Luminex
complexity simplified.

support@luminexcorp.com or orders@luminexcorp.com

For Laboratory Use. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence.

©2016-2018 Luminex Corporation. All rights reserved. Luminex, ARIES, and MultiCode are trademarks of Luminex Corporation, registered in the U.S. and other countries. SYNCT is a trademark of Luminex Corporation.

www.luminexcorp.com

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europe@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@luminexcorp.com

JAPAN

+81.3.5545.7440

info@luminexcorp.com