

INTRODUCING

 GLIDESCOPE® CORE™



SEE MORE.  
DO MORE.

INTRODUCING  
GLIDESCOPE CORE, THE MOST  
COMPREHENSIVE AND FLEXIBLE  
AIRWAY VISUALIZATION SYSTEM  
FOR VIDEO LARYNGOSCOPY AND  
MULTIMODAL AIRWAY PROCEDURES.


GlideScope® Core™ is an all-in-one system offering immediate access to the tools you need to visualize the airway. Designed around a 10" high-definition, touchscreen monitor and comprehensive workstation, GlideScope Core delivers elevated visibility and improved workflow. Discover how GlideScope Core gives you the power to **see more—and do more.**



# ALL NEW GLIDESCOPE CORE

## SEE MORE


GlideScope® Core™ is the first and only airway visualization system to offer live picture-in-picture imaging, giving you the power to view the airway via a GlideScope video laryngoscope and BFlex™, a single-use flexible bronchoscope\*—*simultaneously*.



## GLIDESCOPE CORE PREMIUM FEATURES

- 10" high-definition touchscreen monitor
- Dual connection ports for two video inputs
- Picture-in-picture views—*simultaneously*
- Still image and video capture
- Patient notes annotation feature
- Video playback and image gallery
- SpO2 and pulse rate reading on screen, and on captured videos
- 180 degree image rotation for face-to-face intubations and training
- HDMI out port for external video display

\*GlideScope BFlex, a single-use flexible bronchoscope, is pending premarket notification (510(k)) clearance. Verathon, Inc. cannot be certain that the FDA will or will not grant (510(k)) clearance.



10" TOUCHSCREEN

ADJUSTABLE ARM

CABLE ORGANIZER

BFLEX STORAGE HOOK

PREP TRAY

STERILE BIN

DIRTY BIN

360° ROTATING WHEELS

**GLIDESCOPE CORE WORKS WITH THESE PRODUCTS**

- Spectrum™ single-use blades
- Titanium reusable blades
- BFlex™ single-use flexible bronchoscope\*
- Video Baton 2.0 Large with single-use stats

\*GlideScope BFlex, a single-use flexible bronchoscope, is pending premarket notification (510(k)) clearance. Verathon, Inc. cannot be certain that the FDA will or will not grant (510(k)) clearance.

# DO MORE


## WORKSTATION

The ideal workstation is one that lets you focus more on your patient and less on your equipment. The GlideScope® Core™ system adapts to your space and offers the flexibility to fully focus on your patient. The tools are all neatly organized and close at hand, while the adjustable arm ensures the monitor is always in clear view.


- Mounted on an adjustable arm, the GlideScope Core monitor can be positioned over your patient for optimal viewing in your direct line of sight
- Convenient prep tray ensures easy access to your airway tools
- Sterile bin organizes your GlideScope blades, excess cables and airway management tools in one place to cover a wide range of patients from neonates to large adults
- Dirty bin is designed to support infection control initiatives and house all used equipment for easy disposal
- The organizer attachment keeps your cables ready by the monitor so you can connect your blade easily and start the procedure quickly

## TOOLS

GlideScope Core is the first airway visualization system with dual connection ports and simultaneous live picture-in-picture imaging for GlideScope video laryngoscopes and BFlex™ bronchoscopes\*. No more having to use multiple systems, and vendors, for video laryngoscopy and multimodal airway procedures. With GlideScope Core you can standardize with one system to navigate routine, difficult and emergent airways.



ADJUSTABLE ARM POSITIONS  
EASILY OVER PATIENT



DUAL CONNECTION PORTS WITH  
QUICKCONNECT™ TECHNOLOGY

WHEN IT COMES TO YOUR  
AIRWAY VISUALIZATION NEEDS,  
THE MORE TOOLS YOU HAVE AT  
HAND THE BETTER.

\*GlideScope BFlex, a single-use flexible bronchoscope, is pending premarket notification (510(k)) clearance. Verathon, Inc. cannot be certain that the FDA will or will not grant (510(k)) clearance.

# COMPATIBILITY

GlideScope® Core™ is compatible with the full line of GlideScope visualization tools:

- Spectrum™ single-use blades
- BFlex™ single-use flexible bronchoscope\*
- Titanium reusable blades
- GVL single-use stats



SPECTRUM SINGLE-USE BLADES



BFLEX SINGLE-USE FLEXIBLE BRONCHOSCOPE\*



TITANIUM REUSABLE BLADES



GVL SINGLE-USE STATS



\*GlideScope BFlex, a single-use flexible bronchoscope, is pending premarket notification (510(k)) clearance. Verathon, Inc. cannot be certain that the FDA will or will not grant (510(k)) clearance.

Not all Verathon Inc. products shown are available for commercial sale in all countries.

[Learn more at GlideScopeCore.com](http://GlideScopeCore.com)



[verathon.com](http://verathon.com)

Verathon Inc.  
20001 North Creek Parkway  
Bothell, WA 98011  
USA  
Tel: +1 800 331 2313  
(USA & Canada only)  
Tel: +1 425 867 1348  
Fax: +1 425 883 2896

EC REP

Verathon Medical (Europe) B.V.  
Willem Fenengastraat 13  
1096 BL Amsterdam  
The Netherlands  
Tel: +31 (0) 20 210 30 91  
Fax: +31 (0) 20 210 30 92

Verathon Medical (Australia) Pty Limited  
Unit 9, 39 Herbert Street  
St Leonards NSW 2065  
Australia  
Tel: +1 800 613 603 (Within Australia)  
Fax: +1 800 657 970  
Tel: +61 2 9431 2000 (International)  
Fax: +61 2 9475 1201