

enteral

saf system

SAF System NG Range



Fully NPSA
Compliant



DEHP
Free



Latex
Free



Support
Safe Practice

NEW!

NUTRICARE
LONG TERM
PU NG TUBE

NEW!

NG TUBE
INSERTION
SAFETY
PACK

NEW!

CAREFEED
INFANT
FEEDING
TUBE

CAREFEED
SHORT TERM
NG TUBE

BRIDLE
NG TUBE
RETAINING
SYSTEM

NEW!

ENTERAL
SAF-
EXTENSIONS

NUTRICARE
LONG TERM
PU NG TUBE

NEW!

**FULLY
NPSA
COMPLIANT**


- Fully NPSA Compliant
- Unique Visi-Tip
- **Award winning
NG Safety Insertion
material included**
- Safety Guidewire
- Fully Radiopaque
- Numbered
Every CM

enteral


Safety Guidewire



UNIQUE
VISI-TIP



**Award
winning**
NG Safety
Insertion
material
included



Purple cap
for long
term use.

General

INITIAL
PLACEMENT
NG TUBE

SAFETY
PACK

Placing NG tube

Informed Consent or best
of decision made

Equipment
Patient

X Measurement
on Diagram (1)

and follow traffic light

Guidance if NG tubes
or gastric guidewire not needed

Confirmation, Refer to local policy

ment and on-going checks
recorded in patients notes
policy

Consent Report

RRR001 states;

not be introduced down

gastric placement has

the tube before gastric

men confirmed

stylets should NOT

be placed

gastric placement

ed

Step 2

✓ NE

CHECK

✓ pH value obtained

to feed as per

IT IS

TO F

NO ASP

- No colour

- Change pat

- Perform p

- Flush w

- Offer drink if patient

- Wash all bu

- Advance

Problem solving

RE-ASP

EACH



Colour Coded Tether

FULLY
NPSA
COMPLIANT

USE UP TO
90
DAYS



FULLY
RADIO
OPAQUE

UNIQUE
VISI-TIP

Order Codes - Nutricare-Saf

Colours Codes	Length	Direct Code	NHS Supply Chain Code
6FR	50cm	NC06.050SAF	FWM1876
	80cm	NC06.080SAF	FWM1878
8FR	50cm	NC08.050SAF	FWM1879
	80cm	NC08.080SAF	FWM1880
	92cm	NC08.092SAF	FWM1881
	120cm	NC08.120SAF	FWM1882
10FR	92cm	NC10.092SAF	FWM1883
12FR	92cm	NC12.092SAF	FWM1884

SAF-Extensions

Device Connections	75 cm	150 cm	Direct Code	NHS Supply Chain Code
	●		ESCF.075-SAF	FWM1885
		●	ESCF.150-SAF	FWM1886
		●	EYCF.075-SAF	FWM1887
	●		EYCF.150-SAF	FWM1888

Coming soon quarter four 2013

Reduce the chance of NG Tube Never Events

Increase the adoption of pH testing as 1st line confirmation of NG Placement.

Reduce the incidence and cost of inappropriate X-Rays.

Supports staff inserting NG Tubes to follow NPSA guidelines.

NEW!

NG TUBE INSERTION SAFETY PACK



Dramatically
**Reduce
Pullouts**
with the
AMT Bridle™

BRIDLE
NG TUBE
RETAINING
SYSTEM

Avoid **Pullout Costs**

When NG tubes don't stay in place, costs can mount quickly. Published reports suggest pull outs occur in approximately 40% of cases. Costs of NG tube reinsertion include:

- Cost of clinicians' time to replace NG tube
- Cost of new NG tubes
- Cost to confirm placement via x-ray or fluoroscopy
- **Cost of unnecessary PEG placement**



5-12Fr



14-18Fr

For videos and
support materials
visit our website

www.gbukenteral.com

Carefeed® Infant Feeding Tubes

NEW!

CAREFEED INFANT FEEDING TUBE

Radiopaque
Stripe

Tubing
Features

Clear Tubing

New easy
to open cap
with colour
coded tether.

Numbered every
centimetre
for measured
accurate
placement.

Carefeed® Infant Feeding Tube Sizes

NG Tube Size Length (cm)	4 Fr	5 Fr	6 Fr	8 Fr	10 Fr	12 Fr	14 Fr
48 cm	●	●	●	●			
75 cm	●	●	●	●	●		
92 cm							
125 cm		●	●				

Order Codes - Carefeed® Infant Feeding Tube

Colours Codes	Length	Direct Code	NHS Code
4FR	48cm	CF01.04.048SAF-i	FTA245
	75cm	CF01.04.075SAF-i	FTA246
5FR	48cm	CF01.05.048SAF-i	FTA247
	75cm	CF01.05.075SAF-i	FTA248
	125cm	CF01.05.125SAF-i	FTA249
6FR	48cm	CF01.06.048SAF-i	FTA251
	75cm	CF01.06.075SAF-i	FTA252
	125cm	CF01.06.125SAF-i	FTA253
8FR	48cm	CF01.08.048SAF-i	FTA254
	75cm	CF01.08.075SAF-i	FTA255
10FR	75cm	CF01.10.075SAF-i	FTA256

Radiopaque line easy to check position using X-ray.

Soft, round tip for atraumatic introduction. Countersunk side eyes, for optimal feed delivery.

How to order?

info@gbukenteral.com



Carefeed® Feeding & Drainage Tubes

**CAREFEED
FEEDING &
DRAINAGE
TUBE**

**Tubing
Features**
Fully Radiopaque Tubing

Fully radiopaque easy to check position using X-ray.

Numbered every centimetre for measured accurate placement.



New easy to open cap with colour coded tether.

Carefeed® Feeding & Drainage Tube Sizes

NG Tube Size Length (cm)	4 Fr	5 Fr	6 Fr	8 Fr	10 Fr	12 Fr	14 Fr
48 cm	●	●	●	●			
75 cm	●	●	●	●	●	●	
92 cm				●	●	●	●
125 cm		●	●	●			

Order Codes - Carefeed® Feeding & Drainage Tube Sizes

Colours Codes	Length	Direct Code	NHS Code
4FR	48cm	CF01.04.048SAF	FWM1254
	75cm	CF01.04.075SAF	FWM1255
5FR	48cm	CF01.05.048SAF	FWM1256
	75cm	CF01.05.075SAF	FWM1257
	125cm	CF01.05.125SAF	FWM1258
6FR	48cm	CF01.06.048SAF	FWM1259
	75cm	CF01.06.075SAF	FWM1260
	125cm	CF01.06.125SAF	FWM1261
8FR	48cm	CF01.08.048SAF	FWM1262
	75cm	CF01.08.075SAF	FWM1263
	92cm	CF01.08.092SAF	FWM1615
	125cm	CF01.08.125SAF	FWM1264
10FR	75cm	CF01.10.075SAF	FWM1265
	92cm	CF01.10.092SAF	FWM1616
12FR	75cm	CF01.12.075SAF	FWM1266
	92cm	CF01.12.092SAF	FWM1617
14FR	92cm	CF01.14.092SAF	FTA257

How to order?



+44 (0)1757 282 945

gbukenteral.com
T: +44 (0)1757 282 945

FULLY
NPSA
COMPLIANT
AND DEHP
FREE

enteral

GBUK Enteral Ltd

Blackwood Hall Business Park,
North Duffield, Selby,
North Yorkshire,
YO8 5DD

T: +44 (0)1757 282 945

F: +44 (0)1757 600 545

info@gbukenteral.com

www.gbukenteral.com



Bridle is a registered Trademark
of Applied Medical Technology Inc.
AMT Bridle is manufactured by
Applied Medical Technology Inc.