

ARIES® SARS-CoV-2 Assay (RUO)




The ARIES® SARS-CoV-2 Assay is a qualitative, multiplex, real-time PCR-based test for the detection of SARS-CoV-2 nucleic acid in nasopharyngeal swab (NPS) samples.

Features and Benefits:

- **Full Automation:** A sample-to-answer test that enables targeted SARS CoV-2 detection.
- **Fast Turnaround Time:** Minimal hands-on time and an automated workflow deliver results in about 2 hours.
- **Economical Design:** Just insert your sample into the cartridge—extraction and amplification reagents are included.
- **Robust Detection:** Exonuclease-sensitive probes for the ORF1ab and N viral genes provide exceptional specificity.

Workflow



Ordering Information

Product Name	Part Number
ARIES® SARS-CoV-2 Assay* Protocol File	50-10048 Provided Separately*
ARIES® Two Module System (IVD)	ARIES-M12V1-IVD
ARIES® M1 System (IVD)	ARIES-M6V1-IVD
SYNCT™ Software	CN-SW47

* For Research Use Only.

+ Contact your local Luminex representative or email us at support@luminexcorp.com.

Luminex®
complexity simplified.

orders@luminexcorp.com or support@luminexcorp.com

For Research Use Only. Not for use in diagnostic procedures. Products are region specific and may not be approved in some countries/regions. Please contact Luminex at support@luminexcorp.com to obtain the appropriate product information for your country of residence.

©2020 Luminex Corporation. All rights reserved. ARIES is a trademark of Luminex Corporation, registered in the US and other countries. SYNCT is a trademark of Luminex Corporation.

luminexcorp.com

HEADQUARTERS

UNITED STATES

+1.512.219.8020

info@luminexcorp.com

EUROPE

+31.73.800.1900

europa@luminexcorp.com

CANADA

+1.416.593.4323

info@luminexcorp.com

CHINA

+86.21.8036.9888

info@luminexcorp.com

JAPAN

+81.3.5545.7440

info@luminexcorp.com

SS247721