

Thrombocyte Shaker

Keeps your thrombocytes in motion



Thrombocyte Shaker 96



Thrombocyte Shaker 60

Storage of thrombocytes – standstill or motion

- Upper section moves over 30mm at ~60 strokes at 50Hz / ~72 stokes at 60Hz per minute
- Every shelf has a capacity of 6 bags (for models 24, 48, 60)
- Every shelf has a capacity of 12 bags (for model 96)
- Built in drawer stopper prevents accidental falling out
- Designed in stainless steel
- Each drawer can be operated independently from other drawers
- Designed for silent and reliable long time use
- No-motion alarm
 - Acoustic warn signal
 - D-sub connector for relaying alarms



Thrombocyte Shaker 48



Thrombocyte Shaker 24

Technical data	Thrombocyte Shaker 24	Thrombocyte Shaker 48	Thrombocyte Shaker 60	Thrombocyte Shaker 96
Capacity (bags):	24	48	60	96
Dimensions (cm):	38,1 x 44,8 x 25,9	38,1 x 44,8 x 38,8	38,1 x 44,8 x 45,3	33,8 x 87,7 x 39,3
Weight (kg):	20	25,5	28,25	37,2
Voltage	110-120V\220V-240V and 50HZ/60HZ			
Compliance data				
Classification	Medical device class I			
Directive and Standards applied	MDD 93/42/EEC EN 60601-1; EN 60601-1-2 ; RoHS2			
Classification (Protection against electric shock)	Class I, Type B			