


Burkholderia cepacia complex Now Available in EZ-Accu Shot™ and Epower™

Protect your customers from this notorious environmental contaminant



USP <60> Recommended strains in EZ-Accu Shot™ for Growth Promotion Testing

- Delivers 10-100 CFU per 0.1 ml
- Instant dissolve, no pre-incubation
- Traceability to reference culture ensures authenticity
- Easy & economical refrigerated storage
- Quantitative microorganism preparations delivering 10-100 CFU per 0.1 ml inoculum

Catalog Number and Strain

0488A	<i>Burkholderia cepacia</i> derived from ATCC® 25416™*
01269A	<i>Burkholderia cenocepacia</i> derived from ATCC® BAA-245™*
01270EGPT	<i>Burkholderia multivorans</i> derived from ATCC® BAA-247™
0484A	<i>Pseudomonas aeruginosa</i> derived from ATCC® 9027™*
0485A	<i>Staphylococcus aureus subsp. aureus</i> derived from ATCC® 6538™*

USP <60> Recommended strains in Epower™ for Suitability Testing

- Delivers 1.0-9.9E+07 CFU per pellet
- Ready-to-use format to save you time and money
- Traceability to reference culture ensures authenticity
- Easy & economical refrigerated storage
- Quantitative microorganism preparations delivering 1.0-9.9E+07 CFU per pellet

Catalog Number and Strain

0488E7	<i>Burkholderia cepacia</i> derived from ATCC® 25416™*
00000	<i>Burkholderia cenocepacia</i> derived from ATCC® BAA-245™
01270E7	<i>Burkholderia multivorans</i> derived from ATCC® BAA-247™*
0484E7	<i>Pseudomonas aeruginosa</i> derived from ATCC® 9027™*
0485E7	<i>Staphylococcus aureus subsp. aureus</i> derived from ATCC® 6538™*

Required testing as of December 1, 2019.

***Call us at 1.800.599.BUGS(2847) to place your order.**



*Look for the ATCC Licensed Derivative® Emblem for products derived from ATCC® cultures. The ATCC Licensed Derivative Emblem, the ATCC Licenses Derivative Word mark, and the ATCC Catalog Marks are trademarks of ATCC. Microbiologics, Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.