


Understanding The Difference

MINI Classic Balloon Button VS. MiniONE® Balloon Button



MINI Classic Device

MINI Classic Device

- White Feeding Port
- Non-compatible with Mic-Key® Feeding Sets



MINI Classic Feeding Sets (Seen Above)

- Not compatible with Mic-Key® Device
 - Not compatible with MiniONE® Buttons
- The MINI Classic Button users enjoy the comfort of its low profile and soft, flexible bolster.

MiniONE Device

MiniONE Balloon Button

- Glow Green Feeding Port
- Compatible with Mic-Key® Feeding Sets



MiniONE® Balloon Feeding Sets (Seen Above)

- Compatible with MIC-KEY® Device
 - Non-compatible with Mini Classic Buttons
- The MiniONE offers all the advantages of the MINI Classic but with more features and greater compatibility.
- We recommend the MiniONE Button to all new users as well as previous users.

The Details

Example:

MINI Classic Balloon Button
 14F x 2.3cm = **5-1423**

MiniONE® Balloon Button
 14F x 2.3cm = **M1-5-1423**

- "F" = Tube Diameter (French Size)
- "cm" = Stoma Length (Centimeters)
- Order # Example

MINI Classic	MiniONE	
Yes	Yes	Unique "apple" shape to cover more surface area for a better seal and reduced leakage
Yes	Yes	Soft flexible bolsters for more comfort
No	Yes	Includes size specific Introducer/Stiffener for easier placement
No	Yes	Feed Sets are provided sterile
No	Yes	Feeding Port is Mic-Key® Feed Set compatible
Yes	Yes	Includes comprehensive Patient Education Guide to provide guidance and support
No	Yes	Feeding Port accepts a slip tip syringe for easier medication administration
No	Yes	Glow Green Feeding Port & Feed Set Connector