



TECHNOPATH
CLINICAL DIAGNOSTICS

Multichem[®] ID

INFECTIOUS DISEASE QC

POWERED BY



Multichem[®] ID-SeroNeg

Providing Instrument Specific Third-Party
Infectious Disease Quality Controls



Features

- ◆ Ready to use liquid stable controls
- ◆ Direct sampling from control tubes
- ◆ Instrument specific barcodes
- ◆ Technical and scientific support

Specifications

- ◆ 36 month closed tube stability at 2°C to 8°C.
- ◆ 90 day open tube stability at 2°C to 8°C.

Product	Description	Config.	Part Code
Multichem ID-SeroNeg	Negative Control	4 x 5mL	SR100N
Multichem ID-SeroNeg Mini Kit	Negative Control Mini Kit	1 x 5mL	SR100MN

For further information see Instructions For Use.

THE QUALITY CONTROL COMPANY

Multichem[®] ID-SeroNeg

Analyte List

Negative for the following Analytes;

Anti-HIV IgG

Anti-HCV IgG

Anti-HBc IgG

Anti-HTLV-I IgG

Anti-Treponema Pallidum IgG

HBsAg

HIV p24



TECHNOPATH
CLINICAL DIAGNOSTICS

TECHNOPATH CLINICAL DIAGNOSTICS
Technopath Life Sciences Park, Fort Henry, Ballina, Co. Tipperary, V94 FF1P, Ireland.
For Multichem QC Support: qcsupport@technopathcd.com Website: www.technopathcd.com